# EntraFlo™

NUTRITION DELIVERY SYSTEM

**OPERATING INSTRUCTIONS**  
**GUIDE D'UTILISATION**  
**INSTRUCCIONES DE OPERACIÓN**

**REF** ENT199235

## TABLE OF CONTENTS

(Enteral Delivery Pump  
with Dose Limit and Memory)

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

**CAUTION:** Federal (USA) law and similar laws throughout the world restrict this device to sale by or on the order of a licensed healthcare practitioner.

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
**ES** Números de serie a partir del 850000

### EntraFlo™ Feeding Pump

Model #199235 EntraFlo™ Nutrition Delivery Pump with Dose Limit

### EntraFlo™ DELIVERY SETS

ENFIT70552 1000 mL Bag Pump Set with Pre-attached ENFit™ connector

ENFIT70550 Safety Spike Pump Set with Pre-attached ENFit™ connector

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



---

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---

Before using the EntraFlo™ Nutrition Delivery System read this entire manual. It is important that you understand and are able to follow all the instructions in this manual. If you have questions or concerns about any of the instructions, see back cover to contact your Medline Sales Representative or Customer Service.

## INDICATIONS FOR USE

This pump is intended for the delivery of enteral nutrition and hydration at controlled rates to patient's gastrointestinal system.

## INTRODUCTION

### GENERAL

Upon receipt of a new pump, plug it in for a minimum of 12 hours to charge the backup battery. For maximum battery life, charge the battery at least every 30 days, even if the pump is not in use. Whenever the backup battery is used to power the pump, charge for 12 hours as soon as conveniently possible. The backup battery is designed to provide temporary power to the pump if AC outlet power is interrupted. Frequent, extended operation of the pump on the backup battery will shorten battery life, requiring service by an authorized service representative. The EntraFlo™ Nutrition Feeding Pump is an easy-to-use, versatile device for the appropriate administration of enteral feeding solutions.

### KEY FEATURES OF THIS DEVICE INCLUDE:

- **Ease of Use**
- **Versatility**  
Rate adjustable in 1 mL/hr increments from 1 to 295. Dose adjustable in 5 mL increments from 0 - 9995.
- **Lightweight and Portable**  
An internal backup battery can provide power for short-term use. While the pump should normally be operated on AC outlet power, the internal backup battery can provide operation for up to 8 hours at 100 mL/hr if necessary.
- **Ease of Maintenance**  
The fire resistant plastic case and sealed control panel helps prevent fluid leakage and allow cleaning with any mild soap.
- **Alarms**  
Low Battery, Battery On, Dose Complete, and Occlusion/Empty for detecting empty container or occluded line.
- **½ minute Pause Alarm**
- **Low Occlusion Pressure**  
15 p.s.i. occlusion allows for safe use with any enteral feeding tube.
- **Built-in Pole Clamp**
- **Total Volume Delivered Display**
- **Accumulated Volume Display\***
- **Adjustable Alarm Volume**
- **Dose Limit for Intermittent Feedings**  
With 3 selectable alarm modes.\*
- **Memory Retention of Infusion Rate, Dose Limit, Volume Delivered and Accumulated Volume** even after unit is turned off.

Settings are automatically recalled when the pump is turned back on.

- **Serviceability**  
Self-diagnostic mode allows for complete circuitry check with the push of two buttons.

\*with serial numbers greater than 400000

## CAUTIONS

### PUMP

- NOT FOR PARENTERAL USE INTENDED EXCLUSIVELY FOR ENTERAL FEEDING
- Unauthorized modification to any component of the pump may adversely affect performance or safety in unpredictable ways.
- Turning rotor may pinch fingers.
- Plug and power outlet should be clear of obstruction.
- To disconnect pump from AC power pull plug from power outlet.
- Operator should be within a distance that allows visibility of readouts and ability to hear alarms.
- Wireless communications equipment, such as mobile phones, can affect this device and should be kept at least 4 inches (0.1 meter) from the device during operation.
- Do not operate pump where there is a lot of dust, dirt or airborne particles.
- Do not operate pump with flow sensor in direct path of intense light.
- If necessary to dispose of pump, do so in

accordance with healthcare facility policy and local ordinances for disposal of medical electrical equipment.

- This pump is intended for use only with the available range of EntraFlo™ Administration Sets and only if the Administration Set is correctly installed. Use with other Administration Sets could cause delivery errors and put the patient at risk.
- **Caution:** To reduce the risk of electrical shock, do not remove back cover. Refer servicing to qualified service personnel.
- **Notice:** To insure proper battery operation, unit should be charged for 12 hours prior to use. Remove label prior to use.
- **Danger:** Possible explosion hazard if used in the presence of flammable anesthetics.
- The EntraFlo™ Nutrition Feeding Pump is MEDICAL ELECTRICAL EQUIPMENT requiring special precautions regarding EMC and needs to be installed and put into service according to the EMC information in the ACCOMPANYING DOCUMENTS.
- The EntraFlo™ Nutrition Feeding Pump emits Radio Frequencies (RF) that can affect other MEDICAL ELECTRICAL EQUIPMENT located nearby.

### ADMINISTRATION SETS

- Administration Set tubing may present a risk of entanglement or strangulation. Arrange giving set tubing to minimize the risk of patients or others from strangulation or entanglement.

- Dispose of Administration Set according to healthcare facility procedure for treatment of medical waste, or as directed by a health care provider.
- Follow instructions to reduce the risk of free flow.
- Only for enteral nutrition.

## Alarms/Alerts

**High Priority Alarm** - Flashing red displays and alarm sounds.

Indicates an imminent hazard requiring immediate operator attention, which otherwise could result in serious personal injury and/or product damage if proper action is not taken:

- FREE FLOW

**Medium Priority Alarms** - Flashing yellow indicator light and alarm sounds.

Indicates a condition that requires prompt operator attention, which could result in serious personal injury and/or product damage if proper action is not taken:

- OCCLUSION/EMPTY
- LOW BATTERY

**Alerts** - Indicator light is lit and/or alert sounds  
Indicates that a condition exists of which the operator should be aware:

- LOW BATTERY (light)
- BATTERY CHARGING (light)
- **Note:** Battery Charging light on serial numbers starting at 850000
- RATE CHANGE (alert)
- DOSE COMPLETE (light and alert)

- BATTERY (light)
- PAUSE (alert)

### Cancel Alarm/Alert

Pressing the RUN/PAUSE key after the alarm condition has been corrected cancels the alarm/alert and returns pump to pause mode.

**Note:** See Troubleshooting Section for possible alarm/alert causes.

### FREE FLOW ALARM CONFIRMATION TEST

To determine if the Free Flow Alarm has been installed on the pump, complete the following steps to determine if the pump has been upgraded to include the Free Flow Alarm feature:

1. Turn pump on
2. Press and hold down both the ▲ and ▼ buttons at the same time until the pump beeps
3. Release the buttons
  - a. Below is the information that will be shown:
 

RATE Display	VOLUME DELIVERED Display
888	8888
4. Press and hold down both the ▲ and ▼ buttons at the same time until the pump beeps
5. Release the buttons
  - a. Below is the information that will be shown:
 

RATE Display	VOLUME DELIVERED Display
012	3456
6. Press and hold down both the ▲ and ▼ buttons at the same time until the pump beeps

7. Release the buttons
- a. Below is the information that will be shown:
- RATE Display    VOLUME DELIVERED Display  
4r                      405
8. If pump RATE and VOLUME DELIVERED displays match the information shown in step 7.a., the pump has the Free Flow Alarm feature installed
9. Press and hold the OFF key to exit the Free Flow Alarm test mode

## Symbols

### Serial Numbers less than 700000


<b>IPX1</b>	Degree of Protection Provided by Enclosure
	Medical Electrical Equipment 1. Classified with respect to electrical shock, fire and mechanical hazards in accordance with UL60601-1 2. Classified with respect to electrical shock, fire, mechanical and other specified hazards in accordance with CAN/CSA C22.2 No.601.1

### Serial Numbers 850000 or greater

	Class II Equipment
<b>IP21</b>	Degree of Protection Provided by Enclosure
	Follow Instructions For Use
	Alternating Current
	ON (Power)
	OFF (Power)
	Intertek ETL Classified for Medical Electrical Equipment to IEC 60601-1:2005, IEC 60601-1-11: 2010, CSA C22.2#60601-1:2008

### All Serial Numbers

	Type BF Applied Part (Administration Sets)
	Year of Manufacture
	Attention: Consult Accompanying Documents

	Transportation/Storage Temperature Limits
	Transportation/Storage Humidity Limits
	Manufacturer
	Non-Ionizing Electromagnetic Radiation

## Controls

- ① **ON**, Applies power to unit.
- ② **OFF**, Turns power off except to battery charger.
- ③  Increases flow/infusion rate when pressed and held. Increases dose rate after SET DOSE LIMIT Key has been pressed.
- ④  Decreases flow/infusion rate when pressed and held. Decreases dose after SET DOSE LIMIT Key has been pressed.
- ⑤ **RUN/PAUSE**, Alternately starts and stops unit. Corrects alarm condition and restarts unit when pressed twice.
- ⑥ **RESET**, Resets volume delivered display to 0.
- ⑦ **SET DOSE LIMIT**, Allows dose to be set and displays the dose limit.
- ⑧ **INFUSION RATE DISPLAY**, Indicates delivery rate.
- ⑨ **VOLUME DELIVERED DISPLAY**, Constantly

indicates total volume delivered.

**DOSE LIMIT DISPLAY.** Indicates Dose Limit Volume set when SET DOSE LIMIT Key is pressed during operation.

**ACCUMULATED VOLUME DISPLAY,** indicates Accumulated Volume when SET DOSE LIMIT Key is pressed and held for more than five seconds.\*

- ⑩ **Rotary peristaltic pump mechanism**
- ⑪ **Power cord**
- ⑫ **Alarm volume control**
- ⑬ **Pole clamp**
- ⑭ **Directions**
- ⑮ **Cord wrap**
- ⑯ **Pump adapter bracket**
- ⑰ **Flow sensors and drip chamber guide**
- ⑱ **Tubing guide**

\*with serial numbers greater than 400000

Controls 11; 12; 13; 15 are on the back of the pump

Control 14 is on the side panel

## Operating Instructions

(see Diagrams on page iv)

### DIAGRAM 1:

Activate pump; set desired flow rate.

### DIAGRAM 2:

Insert pump set drip chamber. Make sure set is properly seated on pump.

### DIAGRAM 3:

Assemble silicone tubing around rotor.

### DIAGRAM 4:

Insert pump adapter into bracket.

### DIAGRAM 5:

Place tubing into tubing guide.

Silicone Tubing that is not properly secured around the pump's rotor may result in a free flow situation. Always ensure that the pump adapter is correctly seated into the bracket. This is particularly important when using water or clear liquids, as proper alarm function may not occur if the giving set is not secured properly.

1. Plug the power cord into a 120 VAC., 60 Hz outlet. The internal backup battery should normally only be used for short term use when AC outlet power is unavailable.
2. Press ON - alarm will sound briefly. Remember, the pump has memory and will automatically recall the previous parameter settings. **(DIAGRAM 1)**
3. Check Flow/Infusion Rate. If change is desired, press ▲ or ▼ to change rate setting. **(DIAGRAM 1)**
4. Check Dose Limit by pressing SET DOSE LIMIT. If change is desired, press ▲ or ▼ within the next 3 seconds to change setting limit. Zero Dose Limit implies no dose limit monitoring by pump. **(DIAGRAM 1)**
5. Press RESET to clear Total Volume, if desired.
6. Insert drip chamber into drip chamber

brackets. **(DIAGRAM 2)**

7. Open roller clamp. Expel air from tubing. Close clamp.
8. Stretch silicone tubing around rotor. Be sure not to overstretch the silicone tubing. Insert pump adapter into bridge. **(DIAGRAMS 3+4)**
9. Place tubing into tubing guide. **(DIAGRAM 5)**
10. Attach pump set to patient's feeding tube.
11. Open roller clamp.
12. Press RUN/PAUSE to start pump.
13. OFF will turn off the pump at any time (Flow/Infusion Rate, Dose Limit, Volume\* Delivered and Accumulated Volume are retained.)

\*with serial numbers greater than 400000

## How to Modify Pump Settings

### To change rate or dose limit:

1. Press RUN/PAUSE to stop pump.
2. Select new Flow/Infusion Rate by pressing ▲ or ▼. Select new Dose Limit by pressing SET DOSE LIMIT and then press ▲ or ▼ within next 3 seconds. Press RESET to clear Volume Delivered, if desired. Press RUN/PAUSE to restart pump.

### To view Set Dose Limit:

Press SET DOSE LIMIT. Dose Limit will be displayed for a few seconds.

**To clear volume delivered:**

1. Press RUN/PAUSE to stop pump.
2. Press RESET to clear Volume Delivered.
3. Press RUN/PAUSE to restart pump.

**To view or clear accumulated volume:**

1. Press RUN/PAUSE to stop pump.
2. Press and hold SET DOSE LIMIT for 5 seconds. Total VOL and accumulated volume will be displayed for 10 seconds.
3. To clear accumulated volume, press RESET within the 10 seconds.
4. Press RUN/PAUSE to restart pump.

**To view or change dose complete alarm mode:**

Three alarm modes are available when dose is complete:

- Mode 1 - Pump beeps then stops.
- Mode 2 - Pump beeps, stops and then beeps at 5 minute intervals.
- Mode 3 - Pump beeps continuously.

**To select one of the desired modes:**

1. Press ▲ and ▼ for 5 seconds. The Infusion Rate and Volume Delivered Displays should show 888 and 8888 respectively.
2. Press ▲ and ▼ again. The Infusion Rate Display will show 012 and the Volume Delivered Display will show 3456.
3. Press ▲ and ▼ once again. The Infusion Rate Display will show the current alarm mode (1r, 2r or 3r) (Disregard the Volume Delivered Display).
4. Select desired mode by pressing ▲ or ▼.
5. Press RUN/PAUSE to save new alarm mode.  
**NOTE:** If RUN/PAUSE is not pressed, new

alarm mode will not be saved.

6. Turn pump off.
7. Pump is now set to operate in the new alarm mode.

## Self Test

---

The major functions of the pump may be tested by pressing the ON button and then simultaneously holding the ▲ and ▼ keys for three seconds.

The volume delivered and rate displays should show all eights (8888) and the pump rotor should turn. To test the alarms, press SET DOSE LIMIT, ▲, ▼, RESET or RUN/PAUSE.

When any of these keys are pressed or the flow sensors detect a drop, all alarm lights should light, the alarm should sound and the rotor will halt momentarily.

Any failure indicates a service check should be performed. See Service Section for instructions on reporting a pump in need of service.

## Service

---

EntraFlo™ Nutrition Feeding Pumps are reliable medical electromechanical devices. As with any mechanical device, problems may occur. The following tips will help you correct the minor in-service problems that may arise.

**Note:** Upon receipt of a new pump, plug it in for a minimum of 12 hours to charge the backup battery. For maximum battery life, charge the battery at least every 30 days, even if

the pump is not in use. Whenever the backup battery is used to power the pump, charge for 12 hours as soon as conveniently possible.

<b>Troubleshooting Guide</b>		
<b>Situation</b>	<b>Possible Cause</b>	<b>Solution</b>
BATTERY LIGHT is on	Pump is operating on battery power	Plug cord into AC outlet
LOW BATTERY alarm activated	Battery charge is below performance level	Plug cord into AC outlet
LOW BATTERY light is on	Battery life is fifteen minutes or less	Plug cord into AC outlet To avoid degrading the battery life keep battery charged If low battery persists then return pump for service
OCCLUSION/EMPTY alarm activated	Feeding container is empty	Replace or refill feeding container
	Occlusion is restricting flow	Identify point of occlusion - Pump set · kinked tubing · blocked tubing · closed roller clamp - Feeding tube - Feeding container
	Set drip chamber not properly placed in pump	Assemble properly
	Pump is not in an upright/vertical - position or not clamped properly to the pole	Straighten the pump and clamp properly to the pole
	Pump is on hold for more than 2½ minutes**	Press RUN/PAUSE button
	Bright light is reaching drop sensor	Turn pump so sensor is not in direct light
	Drip chamber walls are coated with formula or condensation	Manipulate chamber to clear drops in sensor path
	Drop sensor lenses in drip chamber guide are dirty	Clean with cotton swab dampened with isopropyl alcohol
	Drip chamber is overfilled	Replace with new Administration Set

Troubleshooting Guide		
Situation	Possible Cause	Solution
Free Flow ("FrE FLO") alarm activated*	Silicone tubing detached from pump rotor; formula is gravity flowing	Assemble silicone tubing properly around pump rotor and secure in pump adapter bracket and tubing guide per diagrams 3, 4, 5 on pages iv
	Drip chamber walls coated with formula or condensation in sensor path	Manipulate chamber to clear drops or condensation in sensor path
	Drip sensor is in direct path of intense light	Turn pump so sensor is not in direct sunlight
	Flow sensors need cleaning	Clean with cotton swab dampened with isopropyl alcohol
DOSE COMPLETE light is on and alert sounds	Pump has delivered preset volume	Press RESET. Press RUN/PAUSE. Pump will operate until dose is completed again
BATTERY CHARGING light is on	AC Power is available and the battery is in charging mode <b>Note:</b> Bat Charging light indicator only available on serial numbers greater than 850000	N/A
Audible Alert sounds when ▲, ▼ or RESET is pressed	When pump is running it is designed to give audible alert when one of these buttons is pressed	N/A

\* For Model #199235 with serial numbers greater than 400000

\*\* Model 199235 only

If the pump serial number is less than or equal to 400000 and the pump RUN/HOLD button is colored RED, refer to section titled Free Flow Alarm Confirmation Test to determine if the Free Flow Alarm has been installed on the pump.

## Cleaning

**Caution:** Always disconnect the pump from the AC outlet before cleaning to avoid electric shock hazard.

Using a cloth or sponge, clean the pump housing and rollers regularly with warm soapy water. (Prolonged exposure to alcohol, household detergents or strong cleaners can result in damage to pump housing.) Clean drip sensors with cotton swabs and isopropyl alcohol.

## Disinfecting

This pump should be disinfected with each change of patient.

If it is necessary to disinfect the entire pump, the following is recommended:

### General disinfectants:

Cidex®  
pHisohex®  
Hibiclens®  
Isopropyl alcohol

### AIDS and hepatitis disinfectants:

10% concentration of 5.25% sodium hypochlorite (household bleach)

### Tuberculosis disinfectants:

70% concentration of isopropyl alcohol

After application of these disinfecting agents, allow to air dry. It should be noted that these recommendations should not be substituted for official procedures that may differ among

institutions. To determine the proper procedure in a particular institution, contact the Infection Control Department.

### DO NOT AUTOCLAVE. DO NOT IMMERSE PUMP IN CLEANING SOLUTIONS.

(Cidex®, pHisohex® and Hibiclens® are trademarks of Johnson and Johnson Corporation, The Mentholatum Company and Regent Medical Limited, respectively.)

## Factory Service

All service work should be performed by authorized personnel only. If the pump does not function properly: Contact your supplier of the EntraFlo™ Feeding Pump to arrange for pump return and replacement. If your pump was supplied directly from Medline contact your Medline Sales Representative or call Customer Service at 1-800-Medline from anywhere in the United States or 1-800-396-6996 from Canada.

## Limited Warranty

Warranty. Company represents and warrants to Buyer that, as of the date of shipment by Company, each and every Product manufactured by Company: (i) is in accordance with the specifications and indications described in the labeling, and (ii) has been properly manufactured, branded, and labeled within the meaning of the US Federal Food, Drug and Cosmetic Act. Company's entire liability and Buyer's exclusive remedy for breach of any warranty shall be, at Company's option, to repair the Product, replace the Product, or provide a credit for the original purchase price. Descriptions or specifications appearing in Company's literature are meant to generally describe the Products and do not constitute any express warranties. The warranties stated herein shall be void and of no effect for any Product failure or deficiency caused by Buyer's shipment, handling, or improper storage; any alteration or modification of the Product; failure to use the Product in accordance with its directions, including use or distribution of the Product after its expiration date; neglect or misuse; or damages occurring after delivery of the Product.

Limitation of Warranty. THE FOREGOING WARRANTIES ARE EXCLUSIVE AND IN LIEU OF ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES WHATSOEVER, AND COMPANY EXPRESSLY DISCLAIMS ALL WARRANTIES OF MERCHANTABILITY, NONINFRINGEMENT, AND FITNESS FOR A PARTICULAR PURPOSE. UNDER NO CIRCUMSTANCES SHALL COMPANY HAVE ANY OBLIGATION OR LIABILITY FOR

ANY LOSS OF PROFIT OR FOR ANY CONSEQUENTIAL, INCIDENTAL, INDIRECT, SPECIAL, PUNITIVE, OR CONTINGENT DAMAGES WHATSOEVER WITH RESPECT TO CLAIMS MADE HEREUNDER OR BY ANY CUSTOMER, CONSUMER OR OTHER USER OF COMPANY PRODUCTS. IN THE EVENT THAT ANY APPLICABLE LAWS IMPOSE WARRANTIES, CONDITIONS, OR OBLIGATIONS THAT CANNOT BE

EXCLUDED OR MODIFIED, THIS PARAGRAPH SHALL APPLY TO THE GREATEST EXTENT ALLOWED BY SUCH LAWS.

## Technical Information

### ELECTROMAGNETIC EMISSIONS

The EntraFlo™ Feeding Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the EntraFlo™ Feeding Pump should ensure that it is used in such an environment.

<b>Guidance And Manufacturer's Declaration-Electromagnetic Emissions For All Equipment And Systems</b>		
<b>Emissions Test</b>	<b>Compliance</b>	<b>Electromagnetic Environment-Guidance</b>
RF Emissions CISPR 11	Group 1	The EntraFlo™ Feeding Pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	
Harmonic Emissions IEC 61000-3-2	Not Applicable	The EntraFlo™ Feeding Pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage Fluctuations/Flicker Emissions IEC 61000-3-3	Not Applicable	

## ELECTROMAGNETIC IMMUNITY

The EntraFlo™ Feeding Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the EntraFlo™ Feeding Pump should ensure that it is used in such an environment.

<b>Guidance And Manufacturer's Declaration-Electromagnetic Immunity For All Equipment And Systems</b>			
<b>Immunity Test</b>	<b>IEC 60601 Test Level</b>	<b>Compliance Level</b>	<b>Electromagnetic Environment Guidance</b>
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV (contact) ± 8 kV (air)	± 6 kV (contact) ± 8 kV (air)	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not Applicable Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Not Applicable Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT (>95% dip in UT ) for 0,5 cycle 40% UT (60% dip in UT ) for 5 cycles 70% UT (30% dip in UT ) for 25 cycles < 5% UT (>95% dip in UT ) for 5 sec	Not Applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EntraFlo™ Feeding Pump requires continued operation during power mains interruptions, it is recommended that the EntraFlo™ Feeding Pump be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<b>NOTE:</b> UT is the a. c. mains voltage prior to application of the test level.			

**Guidance And Manufacturer's Declaration-Electromagnetic Immunity  
For Equipment and Systems That Are Not Life-Supporting**

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conduct RF IEC 61000-4-6	3 V 150 KHz to 80 MHz	Not Applicable	Portable and mobile RF communications equipment should be used no closer to any part of the EntraFlo™ Feeding Pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.17\sqrt{P}$ $d = 0.35\sqrt{P}$ 80 MHz to 800 MHz $d = 0.7\sqrt{P}$ 800 MHz to 2.5 GHz where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). <sup>b</sup> Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol:  
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the EntraFlo™ Feeding Pump is used exceeds the applicable RF compliance level above, the EntraFlo™ Feeding Pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the EntraFlo™ Feeding Pump.

<sup>b</sup> Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3 V/m.

## Separation Distances

### Recommended Separation Distances Between Portable And Mobile RF Communications Equipment And The EntraFlo™ Feeding Pump

The EntraFlo™ Feeding Pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the EntraFlo™ Feeding Pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the EntraFlo™ Feeding Pump as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output Of Transmitter W	Separation Distance According To Frequency Of Transmitter M		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 0.35\sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7\sqrt{P}$
0.01	0.12	0.04	0.07
0.1	0.37	0.11	0.22
1	1.17	0.35	0.70
10	3.69	1.11	2.21
100	11.67	3.50	7.00

For transmitters rated at a maximum output power not listed above the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Specifications

### TYPE OF DEVICE:

Volumetric enteral feeding pump

### MECHANISM:

Rotary Peristaltic

**FLOW/INFUSION RATE:**

1-295 mL/hr in 1 mL increments

**DOSE LIMIT:**

0-9995 mL in 5 mL increments

**ACCURACY:**

±10% of selected flow rate

**OCCLUSSION PRESSURE:**

Does not exceed 15 p.s.i.

**BATTERY:**

- Sealed lead acid
- Rechargeable
- Automatically recharges when connected to an AC outlet
- Recharge time is 12 hours from complete discharge to full recharge
- Operating time is 8 hours at 100 mL/hr (Charge unit fully every 30 days)

The backup battery is designed to provide temporary power to the pump if AC outlet power is interrupted. Frequent, extended operation of the pump on the backup battery will shorten battery life.

**DIMENSIONS:**

7.4" High x 9.5" Wide x 5.3" Deep

**WEIGHT:**

5.7 Pounds

**CASE MATERIAL:**

Fire Resistant Plastic

**POWER REQUIREMENT:**

120V ~, 187mA, 60Hz

**Serial Number less than 700000:**

Protected Against Vertically Falling Water Drops	IPX1
Type of Protection Against Electrical Shock	Class I

**Serial Numbers 850000 or greater:**

Type of Protection Against Electrical Shock	Class II
Protection Against Access to Hazardous Parts with a Finger and Vertically Falling Water Drops	IP21

**All Serial Numbers:**

Expected Battery Life	2 years (200 full charges)
Expected Life of Device	5 years
When integrity of external Mains is in question - use battery power only	
Mode of Operation	Continuous
Operating Temperature	5°C - 40°C (41°F - 104°F) 95% R.H. non-condensing
Storage Temperature	0°C - 50°C (32°F - 122°F) 10% - 95% R.H. non-condensing
Sound Pressure Level	Does not exceed 70dB(A)
Degree Of Protection Against Electrical Shock	TYPE BF (Administration Sets are non-conductive)
Operating Barometric Pressure	700 hPa (10.2 p.s.i.) - 1060 hPa (15.4 p.s.i.) Pump components are not affected by pressure conditions.

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©2016 Manufactured for: Medline Industries, Inc.
EntraFlo is a trademark and Medline is a registered
trademark of Medline Industries Inc. ENFit is a
trademark of GEDSA used with their permission. To