



E Cylinder Oxygen Supply Time Guide

Pressure Gauge	1	2	3	4	5
2000 psi	8 hours	4 hours	2 hours, 30 min	2 hours	1 hour, 30 min
1500 psi	6 hours, 30 min	3 hours	2 hours	1 hour, 30 min	1 hour
1000 psi	4 hours	2 hours	1 hour, 15 min	1 hour	30 min
500 psi	2 hours	1 hour	25 min	15 min	5 min

M6 Cylinder Oxygen Supply Time Guide

	1.5	2	2.5	3	4.5
1	1.5	2	2.5	3	4.5
2.7 hours	1.8 hours	1.4 hours	1.1 hours	0.9 hours	0.6 hours

M60 Cylinder Oxygen Supply Time Guide

	1.5	2	2.5	3	5
1	1.5	2	2.5	3	5
28 hours	21 hours	14 hours	12.25 hours	10.5 hours	5.8 hours