

## Rollator with 6" Wheels



- Comfortable padded seat
- Pouch under seat for added privacy and security of personal items
- Removable, hinged padded backrest can be folded up or down
- 6" Black non-marring casters with soft grip tires are ideal for indoor and outdoor use
- Easy-to-use deluxe loop locks
- Brakes with serrated edges provide firm hold

| PRODUCTS | Description               |
|----------|---------------------------|
| R800BL   | Blue, E0143, E0156, 1/ea  |
| R800RD   | Red, E0143, E0156, 1/ea   |
| R800BK   | Black, E0143, E0156, 1/ea |

| SPECIFICATIONS         | Description                 |
|------------------------|-----------------------------|
| Handle Height          | 32" - 37"                   |
| Unit Dimensions        | 24.5" (W) x 27.5" (L)       |
| Seat Dimension         | 14" (W) x 12" (D) x 20" (H) |
| Weight                 | 16.7 lbs.                   |
| Carton Shipping Weight | 17 lbs.                     |
| Weight Limit           | 300 lbs.                    |
| Warranty               | Limited Lifetime            |

| PARTS        | Description                     |
|--------------|---------------------------------|
| 950340089901 | Handgrip- Lt., 1/ea             |
| 950340079901 | Handgrip- Rt., 1/ea             |
| 9502F80005   | Handlebar Assembly, Left, 1/ea  |
| 9502F80006   | Handlebar Assembly, Right, 1/ea |
| 950300169901 | Anti-rattler, 1/ea              |
| 9501S72602   | Handgrip, 1/pr                  |
| 9501S7261L   | Wheel Lock Left, 1/ea           |
| 9501S7261R   | Wheel Lock Right, 1/ea          |
| 9500002211   | Axle Bearing, 1/ea              |
| 9500001232   | a. Cable Adjuster, 1/ea         |
| 9502M72610   | Fork Assy. 2M Serial #, 1/ea    |
| 9500004101   | c. Cable Fastener, 1/ea         |
| 9501S72600   | Wheel, 1/ea                     |
| 9505W1026120 | Fits in Front of Frame, 1/ea    |
| STDS1060     | Fork Stem Bearing, 1/ea         |
| 9502M72601   | Fork/Caster Assembly, 1/ea      |

| PARTS        | Description                  |
|--------------|------------------------------|
| 9502F80003   | Seat & Brackets, 1/ea        |
| 9502F80004   | a. Seat & Brackets, 1/ea     |
| 9502F80002   | Backrest, 1/ea               |
| 9501S72603   | Backrest Pivot Bracket, 1/pr |
| 9502F80007   | Hand Brake & Cable, 1/ea     |
| 10208C       | Cable Clip, 1/ea             |
| 950330149901 | Axle Cap, 1/ea               |
| 9500001232   | Bushing, 1/ea                |
| 9502F80001   | Folding Hinge Assembly, 1/ea |
| 9502M72610   | Fork Assembly                |
| 9501S72616   | Spring, 1/ea                 |