

## Deluxe Folding Walker, Two Button with 5" Wheels



INSET A

- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter aluminum construction ensures maximum strength while remaining lightweight
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces (Figure A)
- Vinyl-contoured hand grip
- Comes with wheels (Replacement wheels item # 10109)
- The 10210-4ASM Walker comes fully assembled out of the box

### PRODUCTS

PRODUCTS	Description
10210-4	Adult, 4/cs
10211-4	Junior, 4/cs
10210-4ASM	Adult, Assembled, 4/cs
10211-4ASM	Junior, Assembled, 4/cs
10210-1	Adult, 1/cs

### SPECIFICATIONS

SPECIFICATIONS	Description
Width	17" inside hand grip   20" inside back legs
Overall Width	24"
Depth (Opened at Base) (10210   10211   10210-4ASM)	17.5"   16.25"   17.5"
Weight (10210   10211   10210-4ASM)	7.5   7   7.5 lbs.
Carton Shipping Weight (10210   10211   10210-4ASM   10210-1)	30   29   34   9.5 lbs.
Weight Limit	350 lbs.
Warranty	Limited Lifetime