

Heavy Duty 50 PSI Compressor



INSET A



INSET B



INSET C



INSET D

- Newly designed compressor vibration mounts, clamps and lock washers, and wire harnesses improve upon an already durable design
- Ultra-powerful compressor for both humidification and high flow, high-pressure nebulization of viscous medications
- Designed to provide efficient, continuous high-pressure performance
- Lightweight and easy-to-clean

PRODUCTS

Description

18450	Heavy Duty 50 PSI Compressor, 1/cs
-------	------------------------------------

SPECIFICATIONS

Description

Electrical Requirement	115 VAC ~ 60 Hz
Power Consumption	160 W
Sound	~ 60 dBA
Flow Rate	11~41 LPM
Operating Temperature Range	50° - 140° F
Warranty	Limited 1 Year

ACCESSORIES

Description

MASK 001A	Adult Aerosol Mask, 50/cs
MASK 001P	Pediatric Aerosol Mask, 50/cs
MASK 007A	Adult Tracheotomy Mask, 50/cs
MASK 007P	Pediatric Tracheotomy Mask, 50/cs

PARTS

Description

18450-CAP	Filter Cap
18450-CB	Circuit Breaker
18450-FILTER	Intake Filter, 3/pkg
18450-RF	Rear Cabinet Filter
18450-GAUGE	Pressure Gauge
18450-HEADKIT	Compressor Head Rebuild
18450-SWITCH	Power Switch
18450-T	Tee Connector
18450-TUBE	PE Tubing
18450-VALVE	Relief Valve
18450-ADAPTER	Outlet Port Adapter
18450-NVALVE	Pressure Valve Blue Valve
18450-FAN	Cooling Fan