

ALLCARE™ BED



FLOOR HEIGHT

3.6" floor level height (around the width of a smartphone)

Enhances safety, simplifies equipment needs, and its easy on staff

Falls Risk – Floor level (3.6" low height) with optional floor mat to address high risk resident falls. Create an even plane with the safety mat which is intended to help reduce falls, documentation, and potential litigation.

Short Term Rehab – advanced features (cardiac chair positioning, trend/reverse trend), optional scale, optional bed exit to manage higher acuity residents.

Standardization - Improve staff training as well as maintenance with **ONE** bedframe solution across your facility.

With ONE Solution



888-828-7729

MED MIZER®

PROUDLY MADE IN THE U.S.A.

MED-MIZER.COM

Height 3.6"-25" (w/o surface)
 optional 8.6"-30" (w/o surface)

Sleep Surface 35" x 80/84"

Weight Capacity 600 lbs

Options Scale, Bed Extender,
 Safety Mat

ALLCARE™ BED

Falls Prevention



Optional Bi-Fold Safety Mat
 (same height as 6" mattress)

Advanced Positioning



Cardiac Chair



Low Chair



Trendelenburg



Reverse Trendelenburg

MED MIZER®

is at the forefront of innovating medical technology for post acute care. We are focused on partnering with you to help address healthcare challenges in the following key areas—



FALLS PREVENTION



BIARIATRIC CARE



MOBILITY



REHABILITATION



WOUND CARE