



### 3-Position Recliner, Bariatric Extra Wide



FIGURE A



FIGURE D

- 3 positions: upright (Figure A), deep recline and elevated leg rest/footrest (Figure B)
- Retractable lock bar secures the chair in the desired position
- Self-adjusting suspension in the seat prevents bottoming out
- Side panels “pop off” for easy cleaning and maintenance
- Large, blow-molded tray locks in 5 positions and can be stored on the side of the chair
- 5" casters (2 with locks) make it easy to move and secure the recliner
- Wide, comfortable armrests are angled to provide support in any reclined position
- Neutral-colored, sturdy steel frame is easy to clean and maintain
- Moisture barrier on seat prevents seepage
- Comfortable built-in headrest
- Foam and vinyl meet California Technical Bulletin 117

#### PRODUCTS

PRODUCTS	Description
D574EW-BR	Blue Ridge, 1/cs
D574EW-J	Jade, 1/cs
D574EW-R	Rosewood, 1/cs

#### SPECIFICATIONS

SPECIFICATIONS	Description
Overall Dimensions	47.5" (H) x 32" (W) x 41.5" (D)
Width Between Arms	24.5"
Arm Height From Floor	26.5"
Arm Height From Seat	9.5"
Seat Dimensions	21" (H) x 25" (W) x 20" (D)
Carton Shipping Weight	127 lbs.
Weight Limit	500 lbs.
Warranty	Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts

#### PARTS

PARTS	Description
D574EW-1003	1/ea
D574EW-1048	
D574EWBACK-BR	1/bx
D574EWBACK-J	1/ea
D574EWBACK-R	1/bx
D574EWSEAT-BR	1/bx
D574EWSEAT-J	1/ea
D574EWSEAT-R	1/bx
D574EW-1009BR	1/ea
D574EW-1009J	1/ea
D574EW-1009R	1/ea

PARTS	Description
D574EW-1010BR	1/ea
D574EW-1010J	1/ea
D574EW-1010R	1/ea
D574EW-1011	1/ea
D574P-1061	1/ea
D574EW-1035	1/ea
D577P-1012	1/ea
D574EW-1065	1/ea
D574P-1037	1/ea
D574P-1036	1/ea
D574P-1064	1/ea
D574P-1020	1/bx
D574P-1027	1/bx
D574P-1025B	1/bx
D574EW-1017	1/ea
D574EW-1017N	1/ea
D574EW-1070BR	1/ea
D574EW-1070J	1/ea
D574EW-1070R	1/ea
D574EW-1070NL-BR	1/ea
D574EW-1070NR-BR	1/ea
D574EW-1070NL-J	1/ea
D574EW-1070NR-J	1/ea
D574EW-1070NL-R	1/ea
D574EW-1070NR-R	1/ea
D574EW-1013	1/ea
D574EW-1050	1/ea
D574EW-1004	1/ea
6A1208000001	
6A1209000001	
6A1208000018	
D574EW-1003	1/ea
D574P-1048	1/ea
D574P-1048-R	
D574EWBACK-BR	1/bx
D574EWBACK-J	1/ea
D574EWBACK-R	1/bx
D574EWSEAT-BR	1/bx
D574EWSEAT-J	1/ea
D574EWSEAT-R	1/bx
D574EW-1004	1/ea
D574EW-1010BR	1/ea
D574EW-1010J	1/ea
D574EW-1010R	1/ea
D574EW-1011	1/ea
D574P-1061	1/ea
D574EW-1035	1/ea
D577P-1012	1/ea
D574EW-1065	1/ea
D574P-1037	1/ea
D574P-1036	1/ea
D574P-1064	1/ea
D574P-1020	1/bx
D574P-1027	1/bx