

Bariatric Rollator with 8" Wheels



- 8" casters with soft-grip tires are ideal for indoor and outdoor use
- Soft padded oversized seat
- Comes with pouch under seat for added privacy and security of personal items
- Tool-free removable padded backrest for comfort
- Special loop lock made of internal aluminum casting operates easily and ensures safety
- Handles easily adjustable in height with self-threading knob
- Limited Lifetime Warranty

PRODUCTS

	Description
10216BL-1	Blue, E0149+E0156, 1/ea
10216RD-1	Red, E0149+E0156, 1/ea

SPECIFICATIONS

	Description
Handle Height	34" - 39"
Seat to Floor Height	22"
Unit Dimensions	28" (W) x 27.5" (L)
Seat Dimensions	19" (W) x 13" (D)
Weight	23.5 lbs.
Carton Shipping Weight	30.5 lbs.
Weight Limit	400 lbs.
Warranty	Limited Lifetime