



Bariatric Aluminum Folding Walker, Two Button (Wider & Deeper Frame Design) Heavy-Duty Wheels

- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design accommodates individuals up to 500 lbs
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability
- Vinyl-contoured hand grip

PRODUCTS

PRODUCTS	Description
10220-2WW	Adult, 2/cs
10220-1WW	Adult, 1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Width	20.5" inside hand grip 23" inside back legs
Overall Width	27"
Depth (Opened at Base)	23"
Height	32" - 39"
Weight	11 lbs.
Carton Shipping Weight	24 lbs.
Weight Limit	500 lbs.
Warranty	Limited Lifetime

ACCESSORIES

ACCESSORIES	Description
10125HD	Bariatric Walker Tray with Cup Holders, 1/ea