



**Free-Standing Trapeze
with Base (Comparable to Invacare®'s Bariatric Trapeze,
model #BARTRAP)**

- Designed to assist larger individuals in changing positions while in bed and aid in transferring in and out of bed
- Easy-to-assemble tool-free design
- Lightweight design for easy delivery
- Chain is adjustable from edge of headboard to the end of boom by sliding to desired position

PRODUCTS

Description

13049SV	For use by patients up to 650 lbs, 1/cs
---------	---

SPECIFICATIONS

Description

Total Dimensions	72" (H) x 31" (W) x 51" (D)
Inside Width Between Legs	26"
Outside Width Between Legs	31"
Boom Length	37"
Weight	74 lbs.
Carton Shipping Weight	83 lbs.
User Weight Limit	650 lbs.
Warranty	Limited Lifetime