

## Defined Perimeter Mattress Cover



- The Defined Perimeter Mattress Cover creates a raised rail, defined perimeter for enhanced fall prevention, without using patient restraints
- Open in mid-section for safe and easy patient egress and ingress
- Built-in medium density foam bolsters safely and comfortably define edges of bed for patient and prevent falling
- Easy-to-install and secure
- Breathable Nylon makes this a truly universal cover, compatible with both foam and powered air surfaces

### PRODUCTS

	Description
14333	Fits Mattresses 34"-39" (W) x 6"-8" (H) x 80" (L), 1/ea
14333-42	Fits Mattresses 39"-43" (W) x 6"-10" (H) x 80" (L), 1/ea

### SPECIFICATIONS

	Description
Fits Mattresses	34-43" (W) X 80" (L) X 6"-10" (H)
Foam Sections	5½" (W) x 29½" (L) x 4" (H)
Shipping Weight	5 lbs.
Cover Material	Nylon
Foam Density	Medium
Warranty	Limited 1 Year