



Deluxe Bariatric Drop-Arm Commode Assembled

- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode
- Durable, heavy-duty, gray powder-coated steel tubing is sturdy and easy-to-maintain
- Large, durable snap-on seat
- Padded armrests provide comfort and stability
- Pail easily slides in from the front of the commode
- Comes complete with 12 qt commode bucket and cover that removes from front
- Weight Capacity 1000 lbs

PRODUCTS

	Description
11135-1	Commode Oversized DropArm, 1/cs

SPECIFICATIONS

	Description
Width	25.75" between arms 26.25" Outside
Seat Dimensions	23.25" (W) x 18.5" (D)
Height	17.5" - 22"
Weight	21 lbs.
Carton Shipping Weight	26 lbs.
Weight Limit	1000 lbs.
Warranty	Limited Lifetime

ACCESSORIES

	Description
RTL12085	Biodegradable Sanitary Commode Liner, 12/bx, 6 bx/cs

PARTS

	Description
11135A	Arm, 1/ea
11135L-1A	Leg 2A, 1/ea
RTL10389G	Utility ReplaceTips, 1/pr, 8/cs
RTL10389GB	Tips, 1/pr, 10/cs
12002PB	Push Button, 1/ea
111351Z-S	Seat 1Z, 1/ea
11135BRKN	Screws, 1/ea
11135Z-s	Seat, 1/ea
11135B	Backrest, 1/ea
11120B	Commode Pail, 1/ea
11135L	Leg, 1/ea