









## **3-Position Recliner**

Item #	Color	UOM
D574-BR	Blue Ridge	1/cs
D574-J	Jade	1/cs
D574-R	Rosewood	1/cs

- 3-positions: upright (Figure A), deep recline and elevated leg rest/footrest. (Figure B)
- Retractable lock bar secures the chair in the desired position.
- Self-adjusting suspension in the seat prevents bottoming out.
- Side panels "pop off" for easy cleaning and maintenance.
- Large, blow-molded tray locks in 5 positions and can be stored on the side of the chair.
- Specifications

  Overall Dimensions 47" (H) X 25" (W) X 41" (D)

  Width 20" between arms

  Seat to Floor Arm Height 9" 26"

  Seat Dimensions 21" (H) x 19" (W) x 19" (D)

  Carton Shipping Weight 110 lbs.

  Weight Limit 250 lbs.

  Warranty Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts

- 5" casters (2 with locks) make it easy to move and secure the recliner.
- Wide, comfortable armrests are angled to provide support in any reclined position.
- Neutral-colored, sturdy steel frame is easy to clean and maintain.
- Moisture barrier on seat prevents seepage.
- Comfortable built-in headrest. Foam and vinyl meet California Technical Bulletin 117 Standard.