

INVACARE®  
**BAR750**  
BARIATRIC BED

---



---

The **Invacare® BAR750 Bariatric Bed** is a heavy-duty full-electric bed frame capable of supporting up to 750 lb. With the patient in mind, this bed expands from 39" wide to a full 48" and from 80" in length to 88". Four ultra-quiet DC actuators easily reposition the head and foot sections, as well as the bed height. And, for the ultimate in safety, the bed comes standard with a battery backup.

---

## Features

Heavy-duty frame design ensures added strength for client support

Split-spring design eases delivery and setup

Two locking casters on head section and two steering lock casters on foot section

Two high-low motors, one on head section and one on foot section

Pendant allows for dynamic client positioning

Quiet, smooth operation

Integrated battery backup, bed ends and half rails standard

Choose from three mattress width options

BARMATT Mattress options meet fire safety regulations 16CFR1633 and 16CFR1632

650 lb. patient weight capacity

750 lb. total weight capacity

## HCPCS Code: E0302



1. Heavy-duty frame
2. Split spring design
3. Locking casters



Model #		
BAR750 Bariatric Bed, with 1 pair half rails		
Bed Packages		
BARPKG750-1633	BAR750/BARMATT39	
BARPKG750-1-1633	BAR750/BARMATT42	
BARPKG750-2-1633	BAR750/BARMATTEXP	
Accessories		
BARMATT39	Bariatric Mattress, 700 lb. capacity	80" L x 39" W x 6.5" D
BARMATT42	Bariatric Mattress, 700 lb. capacity	80" L x 42" W x 6.5" D
BARMATT48	Bariatric Mattress, 700 lb. capacity	80" L x 48" W x 6.5" D
BARMATTEXP	Bariatric Mattress, 700 lb. capacity	80" x 39" - 48" W x 6.5" D (expandable)
BAR750BOLSTER	Bed extension from 80" to 88"	8" x 39" (expandable)

**Invacare Corporation**  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2012

Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols ™ and  
®. All trademarks are  
owned by or licensed  
to Invacare Corporation  
unless otherwise noted.  
Specifications are subject  
to change without  
notification.  
Form No. 07-010  
Rev. 03/12  
111221



Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Regulation Compliance	Limited Warranty
BAR750IVC	15.5"-29.5" H	80-88" L 39", 42", 48" W	750 lb.	465 lb.	56" L 44" W 41" D	UL962 CSA Certified	5 Years Frame/ Structure 3 Years Mechanical/ Electrical

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

